

# Monitor carts

## FUNCTIONAL ERGONOMICS

Monitor carts are specially designed to accommodate medical monitors which often need to be moved. When selecting the appropriate working height, integrated mechanical or electrical height adjustment are provided for extra comfort. Ergonomic push handles facilitate manoeuvring and positioning. Four lockable twin castors ensure high comfort while moving and a secure fixture when static.



uni-cart

vexio-cart

pro-cart

Monitor carts

# MONITOR CARTS uni-cart



## uni-cart monitor cart with mechanical height adjustment

**Equipment cart, 47 U (Total height: 1685 mm), fitted with:**

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 mechanical height adjustment
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 14 kg / 30.8 lbs
- 1 counterweight set
- 2 cable winders, rear
- 1 push handle, rear

**Installation height:** 47 U (1504 mm)

**Load capacity:** 65 kg / 143 lbs

**Finish:** RAL 7035 light grey powder-coated

**Decorative trims:** push handle: see <sup>1)2)</sup>

**Scope of delivery:** assembled

## uni-cart monitor cart with electric height adjustment

**Equipment cart, 35 U (Total height: 1305 mm), fitted with:**

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 electric height adjustment, load capacity up to 20 kg / 44 lbs with power pack and remote control
- 1 monitor holder, universal with VESA 75 / 100 adaptation up to 14 kg / 30.8 lbs
- 1 keyboard drawer with mousepad (left / right)
- 1 EURO socket strip, 3-way with on/off switch
- 2 cable winders, rear
- 1 push handle, rear
- 1 PC / computer holder (small), rear
- 1 connection cable, Schuko, length 5 m
- 1 integrated plug removal barrier

**Installation height:** 35 U (1120 mm)

**Load capacity:** 65 kg / 143 lbs

**Finish:** RAL 7035 light grey powder-coated

**Decorative trims:** push handle see <sup>1)</sup>

**Scope of delivery:** assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	RS.4914.XXX <sup>1)2)</sup>
10-16 kg / 22-30.8 lbs	RS.4915.XXX <sup>1)2)</sup>

Input voltage	Item number
230 V and 50 / 60 Hz	RS.4145.XXX <sup>1)</sup>
110 V and 50 / 60 Hz	RS.4146.XXX <sup>1)</sup>

<sup>1)</sup> Decorative trims: please replace .XXX with .901 = RAL 5021 water blue, .902 = RAL 5013 cobalt blue, .903 = RAL 7035 light grey

<sup>2)</sup> Decorative trims: please replace .XXX with .800 = RAL 7016 anthracite grey

# MONITOR CARTS vexio-cart



## vexio-cart monitor cart with mechanical height adjustment

**Equipment cart, 23 U (Total height: 933 mm), fitted with:**

- 1 basic frame with 4 locable twin castors Ø 100 mm
- 2 cable ducts, at the side
- 4 hygiene strips
- 1 electric-box in the equipment cart base
- 1 mechanical height adjustment with release pedal
- 1 monitor holder, tiltable with VESA 75 / 100 adaptation up to 14 kg / 30.8 lbs and keyboard shelf with mousepad (left/right)
- 1 PC / computer holder, small
- 1 EURO socket strip, 3-way, with external earthing connector
- 2 cable winders, rear

**Installation height:** 23 U (736 mm)  
**Load capacity:** 65 kg / 143 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** **Support column:** anodised  
**cable ducts:** see <sup>1)</sup>  
**Scope of delivery:** assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	NT.5007.XXX <sup>1)</sup>
8-14 kg / 17.6-30.8 lbs	NT.5008.XXX <sup>1)</sup>

## vexio-cart monitor cart with mechanical height adjustment

**Equipment cart, 23 U (Total height: 933 mm), fitted with:**

- 1 basic frame with 4 locable twin castors Ø 100 mm
- 6 hygiene strips
- 1 mechanical height adjustment with release pedal
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 14 kg / 30.8 lbs with locators
- 1 counterweight set
- 2 cable winders, rear
- 1 push handle, rear

**Installation height:** 23 U (736 mm)  
**Load capacity:** 65 kg / 143 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** **Support column:** anodised  
**push handle / locators** see <sup>1)2)</sup>  
**Scope of delivery:** assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	NT.5014.XXX <sup>1)2)</sup>
8-14 kg / 17.6-30.8 lbs	NT.5015.XXX <sup>1)2)</sup>

# MONITOR CARTS **pro-cart**



## pro-cart monitor cart with electric height adjustment

**Equipment cart, 35 U (Total height: 1332 mm), fitted with:**

- 1 basic frame with 4 locable twin castors Ø 100 mm
- 1 electric height adjustment, load capacity up to 20 kg / 44 lbs with power pack and remote control
- 1 shelf with handle, front
- 1 monitor holder, universal, with VESA 75 / 100 adaptation up to 14 kg / 30.8 lbs
- 1 keyboard drawer with mousepad (left / right)
- 1 push handle, rear
- 2 cable winders, rear
- 1 connection cable, Schuko, length 5 m

**Installation height:** 35 U (1120 mm)  
**Load capacity:** 80 kg / 176 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** handle and push handle: see <sup>1)</sup>  
**Scope of delivery:** assembled

## pro-cart monitor cart with mechanical height adjustment

**Equipment cart, 48 U (Total height: 1798 mm), fitted with:**

- 1 basic frame with 4 locable twin castors Ø 125 mm, conductive
- 1 mechanical height adjustment
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 14 kg / 30.8 lbs
- 1 counterweight set
- 1 push handle, rear
- 2 cable winders, rear

**Installation height:** 48 U (1536 mm)  
**Load capacity:** 80 kg / 176 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** push handle: see <sup>1)2)</sup>  
**Scope of delivery:** assembled

Monitor carts

Input voltage	Item number
230 V and 50 / 60 Hz	<b>PT.5012.XXX</b> <sup>1)</sup>
110 V and 50 / 60 Hz	<b>PT.5017.XXX</b> <sup>1)</sup>

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	<b>PT.5016.XXX</b> <sup>1)2)</sup>
10-16 kg / 22-35.2 lbs	<b>PT.5015.XXX</b> <sup>1)2)</sup>

<sup>1)</sup> Decorative trims: please replace .XXX with .901 = RAL 5021 water blue, .902 = RAL 5013 cobalt blue, .903 = RAL 7035 light grey

<sup>2)</sup> Decorative trims: please replace .XXX with .800 = RAL 7016 anthracite grey

# MONITOR CARTS FOR PHILIPS IntelliVue



## uni-cart rollstand for Philips IntelliVue MP5-70 / MX400-800

**Equipment cart, 32 U (Total height: 1265 mm), fitted with:**

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adapter for Philips IntelliVue MP5-70 / MX400-800 up to 14 kg /30.8 lbs
- 1 vertical holder for „external rack“
- 1 storage basket
- 1 push handle, rear
- 1 integrated plug removal barrier

**Installation height:** 32 U (1024 mm)  
**Load capacity:** 50 kg / 110 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** **push handle** and **holder:** Medium Aqua Accent powder-coated  
**Scope of delivery:** assembled, Philips Table Top Mount is not included in the scope of delivery

Item number

RS.4903.801

## vexio-cart rollstand for Philips IntelliVue MP5-70 / MX400-800

**Equipment cart, 25 U (Total height: 1042 mm), fitted with:**

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adapter for Philips IntelliVue MP5-70 / MX400-800 up to 14 kg /30.8 lbs
- 1 storage basket
- 1 push handle, rear
- 6 hygiene strips

**Installation height:** 25 U (800 mm)  
**Load capacity:** 65 kg / 143 lbs  
**Finish:** RAL 7035 light grey powder-coated  
**Decorative trims:** **push handle:** RAL 7035 light grey powder-coated  
**Scope of delivery:** assembled, Philips Table Top Mount is not included in the scope of delivery

Item number

NT.5016.903

# MONITOR CARTS FOR PHILIPS Efficia



## **vexio-cart** rollstand for Philips Efficia CM 100 / 120 / 150

**Equipment cart, 25 U (Total height: 1042 mm),  
fitted with:**

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adapter for Philips Efficia CM100/120/150 up to 14 kg /30.8 lbs
- 1 storage basket
- 1 push handle, rear
- 6 hygiene strips

**Installation height:** 25 U (800 mm)

**Load capacity:** 65 kg / 143 lbs

**Finish:** RAL 7035 light grey powder-coated

**Decorative trims:** **push handle:** RAL 7035 light grey powder-coated

**Scope of delivery:** assembled, Philips Table Top Mount is not included in the scope of delivery

Item number

**NT.5017.903**